

# UniMove Slide Transfer Sheet Instruction & Safety Manual

Please read this manual before operation, for  
maintenance instructions and safe usage



The enabling power  
of applied knowledge

## Description

UniMove Slide Sheets are used to help prevent manual handling injuries and to eliminate the need for lifting. Made from polyester and coated with silicon, these slide sheet have minimal friction of movement and added durability enabling them withstand infection control washes without warping.

The Slide Sheet is available in two sizes:

- 7890 - Slide transfer sheet, small, 1000 x 1450mm
- 7891 - Slide transfer sheet, large, 2000 x 1450mm

## Warranty

UniMove Warranty covers faulty workmanship or defective material for the original consumer purchase for 12 months from the purchase date. In the event of a defect covered by this warranty,

- We will at our option repair or replace the product.
- The item must be returned to your distributor or selling agent with proof of purchase.
- Use and care instructions have been complied with
- No unauthorised repairs or alterations have been made

## Warning

UniMove Slide sheets are slippery, so please perform transfer slowly making sure patient's head does not hit the top of the bed.

- This product has been designed and manufactured with safety and function in mind for general use.
- For those with special requirements, who may find this product unsuitable, please contact your nearest health professional.





## Use and Care

- We recommend that you seek accredited manual handling training to ensure carers correctly utilise UniMove slide sheets when repositioning patients.
- UniMove Slide sheets can be used singly (folded in half), or with one sheet placed on top of another, depending on the length of transfer and the size of the sheet.
- At all times, it is recommended that at least 2 carers are present to perform patient transfer.
- Never leave a patient unattended on a slide sheet.
- Turn patient gently on to side. Gather up approximately half the doubled sheet and place gathered material under patient as close as possible to their back, with open ends facing direction of travel.
- Turn patient back on to slide sheet. Carer should now be able to pull the gathered material out from the opposite side of patient, so sheet is flat, and patient is now centred on the sheet.
- A carer on either side of the sheet can then grab hold of the top layer of the sheet at the open end and pull carefully in the desired direction with a lunging motion, bending at the knee with a straight back, repositioning the patient.
- For transfers up the bed, place open end of folded sheet towards the bed head.
- For sideways transfers, place open end of folded sheet towards patient and place under patient at least hip to shoulder. Pull in direction of open end.
- Slide sheets are slippery, so please perform transfer slowly making sure patient's head does not hit the top of the bed.
- The aid is a Class 1 Medical Device (non-invasive and non-sterile).

## Cleaning and Maintenance

- To clean your UniMove product hand or machine wash up to 100°C.
- Do not tumble dry, air dry only.
- Regular maintenance inspections are necessary
- Damaged items should be withdrawn from use until rectified

## Important Information

The following is made in lieu of all warranties, expressed or implied.

Sellers and manufacturers only obligation shall be to replace such quality of this product provided to be defective. Neither seller nor manufacturer shall be liable for any injury, loss or damage, direct, indirect or consequential arising out of the inability to use this product.

Before using, the user shall determine the suitability of the product for the intended use, and the user will assume all risk and liability whatsoever in connection therewith. No statement or recommendation not contained herein shall have any force or effect unless in an agreement signed by the representatives of the seller and manufacturer. The patient's nursing care should never be neglected, and when UniMove products are in use, continual assessment should be carried out to ensure the products are still appropriate for the patient's condition.

UniMove is renowned for quality products and design. Due to continuous improvement, our products are subject to change without notice.



PO Box 82 Bentley WA 6982  
T 08 9351 7500 F 08 9350 5299  
[info@unicarehealth.com.au](mailto:info@unicarehealth.com.au)  
[www.unicarehealth.com.au](http://www.unicarehealth.com.au)

Authorised Distributor

